

# adapt<sup>®</sup>200 easiSpec

1

## Seating Assessment / Specification Form

### User Weight






### User Height

Min: 1200mm  
Max: 2100mm

### Chair Action

Asynchronous  
Pelvic-Tilt



- ① Seat to Shoulder (Back Height)
- ② Back of Buttock to Back of Knee (Seat Depth)
- ③ Hip Width at widest point (Seat width)
- ④ Floor to Underside of Knee (Seat Height)

Min	Max
370mm	690mm
330mm	410mm
270mm	410mm
360mm	640mm*

\*Optional Brake Loaded Footrest/Step will allow for any minimum lower leg dimension

Client .....

Company .....

Assessor .....

Date .....

# Seating Assessment / Specification Form

2

## BACK REST

Please enter required dimensions

### CHAIR CODE

Width Required

Height Required

200

WIDTH ..... mm

HEIGHT ..... mm

WIDTH 270 - 410 mm

HEIGHT 370-690 mm

If arms are required please enter dimensions From Seat Surface to underside of elbow

### ARM HEIGHT

WIDTH ..... mm

HEIGHT ..... mm

HEIGHT?

## SEAT

Please enter required dimensions

### CHAIR CODE

Width Required

Depth Required

200

SEATW: to.....mm

SEATD: to.....mm

WIDTH 270 - 410 mm

DEPTH 330 - 410 mm

- Use Back of Buttock to Back of Knee measurement and SUBTRACT 30mm to calculate optimum seat depth
- Seat depths below 350mm will be fixed, above 350mm seat slide depth adjustment of +50mm will be included to a maximum depth of 410mm
- Use Hip to Hip measurement and ADD 50mm to calculate the optimum seat width

## SEAT HEIGHT

Please tick required option (approx heights)

### GAS CODE

MAN 9"

MANUAL LIFT

\*GAS5"

STANDARD GAS COLUMN

\*GAS7"

\*GAS9"

SLOCK7

SWIVEL LOCK

SLOCK9

▲ ON ▼ OFF

## BASE

Please tick required option

Standard Black Base. GAS2 fitted as standard

Polished Aluminium Base

Movement Lock (Non-swivel Gas fitted as standard. Please specify gas lift height required).



## ARM OPTIONS

## Level One Adaptations

<input checked="" type="checkbox"/> Please tick required option	
ARM CODE	HAA*    1234*    3D Trigger    4D*    WSARM    MSA†

## ARM PAD ENHANCEMENTS

<input checked="" type="checkbox"/> Please tick required option	
PAD CODE	UPARM    GEL PADS    *ADARM +50mm    *ADARM -50mm

**ADARM** bracket will add or subtract 50mm from arm height (measured from seat surface to top of arm pad). Applies only to arm options marked with asterisk.

Compatible with HAA & 1234

• Standard seat surface to underside of arm Arm Height is 150 - 240mm. If not suitable use **ADARM** option

## HEAD SUPPORT

<input checked="" type="checkbox"/> Please tick required option	
HEAD CODE	HEAD REST    NECK ROLL    NECK CURVE

## SEAT ENHANCEMENT

## Level Two Adaptations

<input checked="" type="checkbox"/> Please tick required option	
ENHANCE CODE	COCCYX CUT OUT    COCCYX ZONE    EXTRA MEMORY FOAM    THIGH SUPPORT    ERGOCORE CELL /CC    /CZ    /MFS    /TS    /ECORE

<input checked="" type="checkbox"/> Please tick required option	
ENHANCE CODE	WAVE SEAT    SPLIT SEAT    SEWN SEAT WAVE    SPLIT    SEWN SEAT

A sewn seat will provide an even thickness of foam and a uniform tension of fabric across the surface of the seat resulting in a superior sit. This is particularly beneficial when specifying extra layers of memory foam in order to achieve a softer sit

Please specify required **WAVE SEAT DEPTH** ie. Back Buttock to Back of Knee L&R -30mm  
 R (when sitting) .....mm    L (when sitting) .....mm

# **(OPTIONAL) BACK ENHANCEMENTS**


## Level two Adaptations (continued)

adapt<sup>®</sup>200 easiSpec

4

Unless stated otherwise, lumbar cell position will be centered at approx 180mm above seat surface as standard

Please tick box option and state the required dimensions if different from those shown



CODE **LUMBAR CELL** **SACRAL CELL**

Measure distance from seat surface to centre of Air Cell and enter as X+



LUMBAR CELL


LOCATE CENTRE OF AIR CELL AT  
X+ .....mm

Measure distance from seat surface to centre of Air Cell and enter as X+




SACRAL CELL

LOCATE CENTRE OF AIR CELL AT  
X+ .....mm




CODE **/3LUM**

Measure distance from seat surface to centre of Air Cell and enter as X+



/3LUM REACTIVE LUMBAR CELL


LOCATE CENTRE OF AIR CELL AT  
X+ .....mm



CODE **/VSR** **/VSL**


Measure distance from seat surface to centre of Air Cell and enter as X+

Measure distance from centre of back rest to centre of Air Cell and enter as Y+/-



/VSL /VSR VERTICAL CELL


LOCATE CENTRE OF AIR CELL AT  
X+ .....mm  
Y+ .....mm  
Y- .....mm



CODE **/THOR**

Measure distance from seat surface to centre of Air Cell and enter as X+


If not using centrally, measure distance from centre of back rest to centre of Air Cell and enter as Y+/-



THORACIC CELL

LOCATE CENTRE OF AIR CELL AT  
X+ .....mm  
And (if not central)  
Y+ .....mm  
Y- .....mm

Sculpted Back Foam and Bespoke Seat Shaping, please provide separate detailed information



CODE **MFB+SMFB** **SCULP** **SEWN BACK**

Stitched Backrest as standard on SMFB

SOFT WHEEL CASTORS

CODE **SWC**

BRAKE UNLOADED CASTORS

CODE **BUC**

BRAKE LOADED CASTORS

CODE **BLC**